


DIGITAL FULLY AUTOMATIC COLORIMETER

METZ – 402 M



SPECIFICATIONS

Range	400 – 700nm.
Display	Dual Digital Display.
Filters Wave length	400nm,430nm,500nm,520nm,540nm, 580nm,600nm,620nm,680nm.
Display Mode	Absorbance & Wavelength.
Detector	Silicon photo diode.
Zero Adjustment	Auto zero by microprocessor.
Accuracy	± 0.01 O.D.
Amplifier	Micro Controller Based.
Volume	1 ml. Minimum.
Power Supply	230V± 10% AC , 50Hz
Battery Backup	6V Ni – Cd Re-Chargeable Battery.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzerlab.com info@metzerlab.com +91 9897000955