

BINOCULAR METALLURGICAL MICROSCOPE

METZ – 57



SPECIFICATIONS

Magnification	50x – 450x.
Observation	Monocular Head inclined at 45° and rotatable 360°.
Stand	A rigid and stable stand.
Mechanical Stage	Horizontal Mechanical stage 120x140 mm. with fine vernier graduation. Designed with convenient co – axial adjustment for specimen manipulation.
Focusing	Coarse and fine focusing by conveniently placed separate knobs 1 div. 0.002mm.
Illumination	Incident bright light through Epi – illuminator 6V – 20W Halogen lamp with Iris diaphragm and filter slot. Continuously variable luminosity.
Nose Piece	Quadruple revolving Nose piece with positive click stops.
Objectives	Achromatic – M5x, M10x & M45x (SL).
Eye pieces	Wide field 10x.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzerlab.com info@metzerlab.com +91 9897000955